



**NITRILE**

**ACCELERATOR FREE**

**FOOD-GRADE**

*Gloves*



## FEATURE

- Using accelerator free - reducing allergy IV
- Powder-Free, Latex-Free, Silicon-Free
- Textured Fingertips, Beaded Cuff
- Excellent tensile strength and tear-resistance
- Compliance for food-grade requirements
- Non-toxic and suitable for food contact
- Diversified application scenarios
- Using low metal in the raw material & suitable to using in food handling

## PRODUCT INFORMATION



### > Specification

<b>Type</b>	Powder-free and non-sterile
<b>Material</b>	100% nitrile (Acrylonitrile-butadiene). No natural rubber latex
<b>Feature</b>	Ambidextrous, beaded cuff and finger textured surface for easy grip
<b>Color</b>	Blue , White, Black, Purple, Orange, Green, Pink, Yellow
<b>Size</b>	Small, Medium, Large, Extra-large
<b>Donning Agent</b>	The gloves are powder-free and use chlorine to smooth contact surfaces

## APPLICATION SETTINGS

Medical grade & food grade, AQL 1.5, used as good protection to avoid cross contamination between foods, as well as between caterers and foods.

Protection against bacteria and fungi, skin irritation and sensitization tested.



## PHYSICAL PROPERTIES

Property	S&S Glove	ASTM D6319
<b>Elongation at break (%)</b>		
• Before Aging	500 Min	500 Min
• After Aging	400 Min	400 Min
<b>Tensile Stength (Mpa)</b>		
• Before Aging	16 Min	14 Min
• After Aging	16 Min	14 Min

  

Property	S&S Glove	ASTM D455-2
<b>Forece at Break (N)</b>		
• Before Aging	7 Min	6 Min
• After Aging	7 Min	6 Min



## PRODUCT



## STANDARD CARTON SIZE



## CERTIFICATIONS



## STANDARDS

- FDA 510K
- CE PPE CAT III TYPE B
- ASTM D6319
- ASTM F1671
- FDA 21CFR 177
- EN 455
- Standards and Specifications for Food Utensils, Containers and Packages 2021
- Japan Notification 370
- EC 1935/20





## **S&S GLOVE CORPORATION**

Office

Lot Vb.20a2, Road 24, Tan Thuan Export Processing Zone,  
District 7, Ho Chi Minh City, Vietnam

Factory

Lot 4, D6 Street, Dat Do I Industrial Park, Phuoc Long Tho Ward,  
Dat Do District, Ba Ria - Vung Tau Province, Vietnam

Telephone

**(+84) 283 8725 999**

Hotline

**(+84) 987 187 178**

Email

**info@ssglove.vn**

Website

**www.ssglove.vn**



Get in touch with us!



**/ssglovecorp**