

NITRILE MOISTURE *Gloves*



FEATURE

Aloe Vera coating in the donning side that soothes and moistures the skin. A good barrier protection against bacteria, viruses, chemicals and dangerous substances. A dependable protection for tougher jobs for laboratory, cleaning, maintenance and general applications. Free from latex protein, eliminating Type I Immediate Hypersensitivity reaction. Custom design enhanced comfort and fit. Finger textured enhanced wet and dry grip with added handling precision. Chlorinated powder free interior promote smooth donning & doffing experience.



PHYSICAL DIMENSIONS

Dimensions		Standards - S&S Glove				
Length (mm)		Min 230				
Finger Length (mm)	Width	Thumb	Index	Mid	Ring	Little
• Small	85 ± 5	53 ± 5	66 ± 5	77 ± 5	70 ± 5	54 ± 5
• Medium	95 ± 5	55 ± 5	72 ± 5	80 ± 5	75 ± 5	56 ± 5
• Large	106 ± 5	58 ± 5	72 ± 5	82 ± 5	74 ± 5	60 ± 5
• Extra-large	115 ± 5	62 ± 5	78 ± 5	88 ± 5	80 ± 5	62 ± 5

PHYSICAL PROPERTIES

Property	S&S Glove	ASTM D6319
Elongation at break (%)		
• Before Aging	500 Min	500 Min
• After Aging	400 Min	400 Min
Tensile Strength (Mpa)		
• Before Aging	14 Min	14 Min
• After Aging	14 Min	14 Min

QUALITY STANDARDS

Inspection Levels and AQLs

Characteristic	Inspection Level	AQL
Physical dimensions	S-2	4.0
Watertight test	G-1	1.5
Physical Properties	S-2	4.0



SPECIFICATION

Material	100% NBR (synthetic rubber)
Powder content	Powder-free
Color	Green
Thickness	4-6mil
Cuff style	Beading cuff
External glove surface	Smooth with micro-textured finish
Internal glove surface	Aloe vera coated
Allergy prevention	Latex free
Sterile	No

PRODUCT



CERTIFICATIONS



EN455
1, 2, 3



STANDARD CARTON SIZE



STANDARDS

- EN 374
- EN 420:2003 + A1:2009
- EN455
- ASTM 6319
- ISO 10993

MANUFACTURING ACCREDITATIONS

- ISO 9001 Quality Management System Certification
- ISO 13485 Medical Device Quality Management System
- ISO 14001 Environmental Management System
- ISO 45001:2018 Occupational health and safety management systems
- ISO 50001 Energy Management
- ISO/IEC 17025 State Accredited Laboratory





S&S GLOVE CORPORATION

Office

Lot Vb.20a2, Road 24, Tan Thuan Export Processing Zone,
District 7, Ho Chi Minh City, Vietnam

Factory

Lot 4, D6 Street, Dat Do I Industrial Park, Phuoc Long Tho Ward,
Dat Do District, Ba Ria - Vung Tau Province, Vietnam

Telephone

(+84) 283 8725 999

Hotline

(+84) 987 187 178

Email

info@ssglove.vn

Website

www.ssglove.vn



Get in touch with us!



/ssglovecorp