



BIODEGRADABLE

NITRILE

Gloves

FEATURE

- 100% Nitrile, Biodegradable, Powder-Free, Surface-Chlorinated
- Beaded cuff ensures easy donning and prevents rolling down
- Finger textured enhance wet and dry grip
- Decrease risk of allergies
- Superior strength with better puncture resistance
- Excellent tactile sensitivity



PHYSICAL DIMENSIONS

2.5 Mil Biodegradable Nitrile Gloves	
Glove Length (mm/inches)	240/9.5 Min
Palm Thickness (mm/mil)	0.06 ± 0.02/2.4 ± 0.8
Finger Thickness (mm/mil)	0.09 ± 0.02/3.5 ± 0.8
Unit Weight (g)	S: 3.5 ± 0.3 M: 3.7 ± 0.3 L: 4.0 ± 0.3 XL: 4.2 ± 0.3
Color	Green

APPLICATION SETTINGS

Different specifications offer perfect reliable strength and comfortable dexterity for anything from medical healthcare and automotive work to food prep, industrial, janitorial or chemical applications, especially use to protect the environment by it's feature.



PRODUCT



CERTIFICATIONS



STANDARD CARTON SIZE



STANDARDS

- EN455
- EN374
- ASTM D6319
- ASTM D5511
- CE PPE CAT III TYPE B
- CE MARKING

MANUFACTURING ACCREDITATIONS

- ISO 9001
- ISO 13485
- ISO 14001
- ISO 45001:2018
- ISO 50001
- ISO/IEC 17025





S&S GLOVE CORPORATION

Office

Lot Vb.20a2, Road 24, Tan Thuan Export Processing Zone,
District 7, Ho Chi Minh City, Vietnam

Factory

Lot 4, D6 Street, Dat Do I Industrial Park, Phuoc Long Tho Ward,
Dat Do District, Ba Ria - Vung Tau Province, Vietnam

Telephone

(+84) 283 8725 999

Hotline

(+84) 987 187 178

Email

info@ssglove.vn

Website

www.ssglove.vn



Get in touch with us!



/ssglovecorp